

# Weekly Update

Oct 25th, 2017


# Gel Testing

- Tried non stick parchment paper...
- Difficult because can't get the correct shape
- Can't tape it either... non stick..

# Iris

- Tested the irises..
- Looked at transmission with it open, closed, 50%, and 20%
- Not sure how to actually calculate percentage.. Calculating area difference?

Transmission % using Irises



# Next step?

- Iris seems consistent throughout the wavelengths
- Still not sure what to do with the “focussing effect” from earlier...
- Try cuvette? Need more understanding of the detector